
ROSSMAX

www.rossmax.nl

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
### EN English

#### Introduction

The Digital Medical Thermometer provides quick, highly accurate reading over the human medical temperature range. The quality of the thermometer has been verified and conforms to the provisions of the EC council directive 93/42/EEC (Medical Device Directive) Annex I essential requirements and applied harmonized standards. EN 12470-3: 2000 / A1: 2009 Clinical thermometers-Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.

Attention: Consult the accompanying documents. Please read this manual carefully before use. Please be sure to keep this manual.

#### Battery Replacement

Place a new 1.5V DC button size battery type SR41 or LR41 or equivalent in the chamber with positive side faced up and negative side faced down.

#### Display on the LCD

M: Last measured temperature.

■ : When the ■ appears on the bottom right corner of the LCD, the battery is exhausted and needs replacement.

#### Troubleshooting

| Error message                 | Problem                                           | Solution                                                                                                     |
|-------------------------------|---------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Err</b>                    | The system is not functioning properly            | Unload the battery, wait for 1 min and power it. If the message reappears, contact the retailer for service. |
| <b>H 1 °C (bi-bi...bi)x10</b> | Temperature taken is higher than 42.9°C (109.9°F) | Take a new temperature via close contact and sufficient rest.                                                |
| <b>Lo °C</b>                  | Temperature taken is lower than 32°C (90.0°F)     | Take a new temperature via close contact and sufficient rest.                                                |

#### Cleaning and Disinfection

1. Clean the unit by wiping it with a dry cloth and disinfect the probe with ethyl alcohol.

2. Don't let the unit contact any chemical thinner.

3. Do not soak the display in water.

#### Caution

■: Do not bend, drop or twist the thermometer. It is not shockproof.  
Do not store the unit under direct sunlight, at a high temperature, in high humidity or dust.  
Do not disassemble. See BATTERY REPLACEMENT to replace battery.  
Avoid strong electromagnetic interference such as microwave ovens and cell phones.  
Keep battery away from children. Harmful is swallowed.  
Battery should not be charged or placed into extreme heat as it may explode.  
Remove battery from the thermometer when not in operation for a long time.  
If the unit is stored in an extreme environment before measurement, please make sure it has been first adapted in the room temperature, so that the unit can be measured normally.

Measurement results are for reference only. Contact your physician if you have or suspect any medical problems. Do not change your medications without the advice of your physician or healthcare professional.

This device may not meet its performance specification if stored or used outside temperature and humidity ranges specified in specifications.

If this device is used according to the operations instruction, periodic recalibration is not required. If you still have questions, please send the device to dealers.

Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.

High, prolonged fever requires medical attention especially for young children. Please contact your physician.

For safe reason, during children's temperature measurement, please keep them from crying, talking and any related dangerous activities.

Usage of the probe cover may result in a 0.1°C(0.2°F) difference from actual temperature.

To clean the probe before and after using the thermometer to ensure an accurate reading and avoid cross contamination approximately, no matter the probe cover has been installed or not.

Do not drink hot or cold fluids, exercise, and smoke or perform other activities prior to a reading. These activities will raise or lower temperature readings when compared to your normal, average temperature.

Please ensure the close contact between the probe and the armpit of the baby or children.

Do not service or maintain device while in use.

#### Specifications

Measurement Range 32.0°C~42.9°C (90.0~109.9°F)

Accuracy ±0.1°C 35.5°C~42.0°C (±0.2°F 95.9°F~107.6°F); ±0.2°C under 35.5°C or over 42.0°C (±0.4°F under 95.9°F or over 107.6°F)

Display liquid crystal display 3 1/2 digits

Battery One 1.5 V DC button size battery (SR41 or LR41)

Power consumption 0.15 milliwatt in measurement mode

Battery Life Approx. 1500 times operation of 1 year with 1-2 measurement per day including standby mode

Alarm: Approx. 10 seconds sound signal when peak temperature reached

Dimension 12.5cm x 2.1cm x 1.2cm (L x W x H)

Weight Approx. 12.8 grams including battery

Environmental for using Temperature: 5°C~40°C (41°F~104°F), Humidity: 15%~95%RH; 700~1060 hPa

Storage and transport Temperature: 5°C~40°C (41°F~104°F), Humidity: 15%~95%RH; 700~1060 hPa

IP Classification IP22, Protection against harmful ingress of water and particulate matter

Safety classification **1** Type BF Applied Part